 Chest measurement ABOVE breasts:	Standard cup size ranges D, DD, E	Plus cup size ranges F, FF and above
	Order Code:	Order Code:
75 - 84 cm	QADV-XS	QADV-XS Plus
84 - 92 cm	QADV-S	QADV-S Plus
92 - 100 cm	QADV-M	QADV-M Plus
100 - 108 cm	QADV-L	QADV-L Plus
108 - 113 cm	QADV-XL	QADV-XL Plus
113 - 122 cm	QADV-XXL	QADV-XXL Plus
122 - 134 cm	QADV-XXXL	QADV-XXXL Plus
134 - 148 cm	QADV-4XL	QADV-4XL Plus
If the cup size is smaller than D we recommend to use the QualiBreath sternum and thorax support.		
Accessories		Order code
Postoperative care kit (QualiBra Advanced + QualiPad + Axilla Pad)		QTKIT + QADV size
QualiPad seatbelt pad (single unit)		QTSEPA

Does not contain natural rubber latex (no NRL).
Supplied single packed with Instructions for Use.

Washing and care: QualiBra Advanced is washable at 86°F. Do not tumble dry.

To order call +39 0125 54 691

or Email: orders@qualiteam.com

or buy online at www.qualiteam.com

Inside:
Surgeon's skill


Outside:
QualiBra Advanced
protection

QUALIBra
advanced


Nurse-designed
sternum and breast support
for women

Why the internal closure of the sternum needs protection from the outside:


Sternal wound infections:




Increased risk for women and obese patients:




Chronic pain and pulmonary complications:



Average cost of sternal wound complications:



Readmission rates in cardiac surgery are among the highest: 16% on average



Prevent complications, advance recovery, and decrease costs of health care.

QUALIBra
advanced

Postoperative sternum and breast support

Decrease sternal wound infections:

QualiBra Advanced complements the internal closure by providing constant lateral stabilization of the sternum, and the patient can use the integrated pockets to reinforce support at peak internal thoracic pressure (e.g. coughing, sneezing, straining).

Vertical bands at each side of the bra cups prevent the pull of gravity and separate the breast tissue to keep the wound free from buildup of moisture and heat, which helps to protect the sternal wound from infection.



Decrease pain and pulmonary complications:

QualiBra Advanced reduces pain as the vertical bands permit easy, painless access to the wound while maintaining breast support.

The sternum and breast support is a constant reminder to the patient that her wound is being protected and it helps to encourage deeper breathing, early mobilization and supports when coughing.



Advance recovery and esthetic wound healing:

QualiBra Advanced allows the patient to have better control of pain and respiration, which helps to promote earlier mobilization and advance recovery.

The vertical bands keep the breast tissue gently away from the sternal wound which aids in more esthetic healing women of all ages will appreciate.



Improve wound care and protection beyond the hospital:

The sternum band allows for quick and straightforward access to the wound in emergencies, and for easy inspection and care without the need to open the bra or uncover the breasts.

Soft, elasticized cotton bra cups with no seams or labels against the skin, and an ultra-light ventilated chest band together with easily adjustable closures, provide exceptional comfort and enhance compliance for usage at home until the sternum and surgical wound heal.

